

Certificate of Analysis

Product code: OODSF

Product name: OOSAFE® SURFACE DISINFECTANT

Batch/LOT No.: 70944

Production date: 2022-05

Expiry date: 2025-05

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear, slightly yellow liquid	Clear, slightly yellow liquid	30/06/2022
Odour / Olfactive	Slight smell of propionic acid	Slight smell of propionic acid	30/06/2022
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	5.5 – 7.0	6.6	30/06/2022
Density, 20 °C, gr/cm ³ / DIN 51757	1.01 – 1.15	1.03	30/06/2022
Quaternary Ammonium Compounds / QUAT test stripe		<7.0%	30/06/2022
Additives		<6.0%	30/06/2022
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.		SPAR-9128-0722	24/07/2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.		S5969-0722SPAR	21/07/2022

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Andrius Daugintis,
Quality Control department
SparMED ApS
Date: 28/07/2022

