

Certificate of Analysis

Product name: Oosafe® Laboratory Surface Disinfectant

Product code: OODSF

Batch/LOT No.: 70954

Production date: 2022-11

Expiry date: 2025-11

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear, slightly yellow liquid	Clear, slightly yellow liquid	N/A
Odor / Olfactive	slight smell of propionic acid	slight smell of propionic acid	N/A
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	5,5 – 7,0	6,5	N/A
Density, 20 °C, gr/cm ³ / Density measuring device	1,01 – 1,15	1,04	N/A
Quaternary Ammonium Compounds / QUAT test stripe		< 7,0 %	N/A
Additives		< 6,0 %	N/A
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR – 1108 – 0123	Passed	13/01/2023
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S6284 – 0123SPAR	Passed	06/01/2023

Quality Control according to ISO 13458:2016

GOosafe with SparMED!


Diana Magnussen


Quality Control Department

SparMED ApS

Date: 16/01/2023

dpmagnussen



 GOosafe in IVF Laboratories	PREPARED BY:	APPROVED BY:
	QUALITY CONTROL MANAGER	MANAGEMENT REPRESENTATIVE
	Signature <i>lvei</i>	Signature <i>Bmbe</i>