




SparMED Aps  
Ryttermarken 2  
3520 Farum  
Denmark



ELI Accession Number: SPAR-9203-0818

Date of completion: 08-10-2018

Lot number: 82504

Reference number: OODW-70

Description of test articles: Oosafe Disinfection Wipes

**Assay system requested by customer:** An incubator was cleaned with the test article. Post cleaning a culture plate was set up and one cell mouse embryos were cultured in the cleaned incubator for 96-hours.

**Control assay method and results:** 15 one cell (B6C3F1 X B6D2F1) embryos were cultured in triplicate micro drops of culture medium in control incubator ELI-182:

15 / 15 (100 %)

1-cell to 2-cell within 24 hr

15 / 15 (100 %)

1-cell to expanded blastocyst within 96 hr

*For a valid assay, Embryotech™ requires at least 70% of one cell control embryos to develop to expanded blastocyst within 96-hours.*

**Test assay method and results:** 21 one cell (B6C3F1 X B6D2F1) embryos were cultured in triplicate micro drops of culture medium while in incubator ELI-158 which was cleaned with the test article:

21 / 21 (100 %)

1-cell to 2-cell within 24 hr

21 / 21 (100 %)

1-cell to expanded blastocyst within 96 hr

**Summary of observations:** All test and control embryos were selected randomly from a common pool of freshly collected embryos. 100 percent of the control embryos developed to the expanded blastocyst stage within 96-hours. 100 percent of the embryos cultured in an incubator previously cleaned with the test article developed to the expanded blastocyst stage within 96-hours.

Signature  
Study Director

Date

Signature  
Quality Reviewer

Date