



SparMED Aps
Ryttermarken 2
3520 Farum
Denmark



140 Hale Street
Haverhill, MA 01830
qc@embryotech.com

ELI Accession Number: SPAR-9774-1118

Date of completion: 11-20-2018

Lot numbers: 84406

Reference numbers: OODW-70 / 460

Description of test articles: Oosafe Disinfection Wipes

Assay system requested by customer: An incubator was cleaned with the test article. Post cleaning a culture plate was set up and one-cell mouse embryos were cultured in the cleaned incubator with the test article for 96-hours.

Control assay method and results: 15 one cell (B6C3F1 X B6D2F1) embryos were cultured in triplicate micro drops of culture medium in control incubator ELI-181:

15 / 15 (100 %)

1-cell to 2-cell within 24 hr

15 / 15 (100 %)

1-cell to expanded blastocyst within 96 hr

For a valid assay, Embryotech™ requires at least 70% of one cell control embryos to develop to expanded blastocyst within 96-hours.

Test assay method and results: 21 one cell (B6C3F1 X B6D2F1) embryos were cultured in triplicate micro drops of culture medium while in incubator ELI-182 cleaned with test article:

21 / 21 (100 %)

1-cell to 2-cell within 24 hr

20 / 21 (95 %)

1-cell to expanded blastocyst within 96 hr

Summary of observations: All test and control embryos were selected randomly from a common pool of freshly collected embryos. 100 percent of the control embryos developed to the expanded blastocyst stage within 96-hours. 95 percent of the embryos cultured in an incubator previously cleaned with the test article developed to the expanded blastocyst stage within 96-hours.

Signature
Study Director

11-20-2018

Date

Signature
Quality Reviewer

11-20-2018

Date