

Certificate of Analysis

Product name: OOSAFE® CLASS IIA MEDICAL DEVICE DISINFECTANT WIPES

Product code: OOMDW-120

Batch/LOT No.: 92303-MD

Production date: 2019-05

Expiry date: 2022-05

| Analysis / Method | Specification / Reference No. | Result | Completion date |
|--|-------------------------------|---------------------------|-----------------|
| Appearance / Visual | Clear liquid inside wipes | Clear liquid inside wipes | NA |
| Odor / Olfactive | Odourless | Odourless | NA |
| pH-Rate (2% active substance in water) / pH meter (DIN EN 1262) | 7,0-9,0 | 7,8 | NA |
| Density, 20 °C, gr/cm ³ / DIN 51757 | 1,001-1,010 | 1,005 | NA |
| Quaternary Ammonium Compounds / QUAT test stripe | 1,9-2,2 % | 2 % | NA |
| Water | | The rest | NA |
| Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours. | SPAR-2056-0619 | Passed | 01-07-2019 |
| Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven. | S3928-0619SPAR | Passed | 02-07-2019 |

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Inga Sileikaite,
Quality Control Department
SparMED ApS




Date: 22-07-2019