



Doc. No.: QC_4.14.1. Rev. Date: 01.04.2019



GOosafe in IVF Laboratories

SparMED ApS

VAT no.: DK 30 89 85 75

Ryttermarken 2, 3520, Farum, Denmark

Phone: +45 39 40 25 03

e-mail: info@sparmed.dk - www.sparmed.dk

Certificate of Analysis

Product name: OOSAFE® CLASS IIA MEDICAL DEVICE DISINFECTANT WIPES

Product code: OOMDW-120

Batch/LOT No.: 93204-MD

Production date: 2019-07

Expiry date: 2022-07

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid inside wipes	Clear liquid inside wipes	NA
Odor / Olfactive	Odorless	Odorless	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7,0-9,0	8,2	NA
Density, 20 °C, gr/cm ³ / DIN 51757	1,001-1,010	1,005	NA
Quaternary Ammonium Compounds / QUAT test stripe	1,9-2,2 %	2,0 %	NA
Water		The rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-2360-0819	Passed	26-08-2019
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S4003-0819SPAR	Passed	23-08-2019

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Inga Sileikaite,
Quality Control Department
SparMED ApS



Date: 30-Aug-2019