GYNEMED Semen Diagnostics

SemenlgA



- Ready-to-use
- Detection of anti-spermantibodies (ASA). Indicates presence of immunglobulin A (IgA)-type antibodies on motile sperm cells. Detection with Anti-IgA-coated red microspheres
- Shelf life is 18 months from time of manufacture if the product are storde at 2-8°C.

Ref. No.	description
ZR11200	1 bottle containing 300 µl latex particles coated with spec. AntilgA antibodies

SemenIgG



- Ready-to-use
- Detection of anti-spermantibodies (ASA). Indicates presence of immunglobulin G (IgG)-type antibodies on motile sperm cells. Detection with Anti-IgG-coated blue microspheres
- Shelf life is 18 months from time of manufacture if the product are storde at 2-8°C.

Ref. No.	description
ZR11400	1 bottle containing 300 µl latex particles coated with spec. AntilgG antibodies

SemenMar



- Ready-to-use
- Detection of anti-spermantibodies (ASA). Indicates presence of IgG and/ or IgA-antibodies on motile sperm cells. Detection with Anti-IgG and Anti-IgA-coated yellow microspheres
- Shelf life is 18 months from time of manufacture if the product are storde at 2-8°C.

Ref. No.	description
ZR11100	1 bottle containing 300 µl latex particles coated with spec. Anti-Mar antibodies

GYNEMED Semen Diagnostics

SemenStain



- Morphology assessment of spermatozoa. Succedanstaining method.
- Differential staining of spermatozoa tissues based on their basophilic, eosinophilic and neutrophilic properties
- Shelf life is 36 months from time of manufacture if the product are storde at 15-25°C.

Ref. No.	description
ZR10050	1 x 50 ml Reagenz 1
	1 x 50 ml Reagenz 2
	1 x 50 ml Reagenz 3
	1 x 50 ml Reagenz 4

SemenVIT



- Vitality assessment of spermatozoa. Dye exclusion method.
- Based on the integrity of the sperm membrane in live spermatozoa
- Shelf life is 24 months from time of manufacture if the product are storde at 2-8°C.

Ref. No.	description
ZR10300	1 x 20 ml Reagenz 1
	1 x 30 ml Reagenz 2

SemenHOS



- Vitality assessment of spermatozoa. Hypoosmotic swelling method
- Based on the semi-permeability of the intact cell membrane and the possibiling of active water transport
- Swelling effect in intactsperm cells can be observed for up to 30 min
- Shelf life is 24 months from time of manufacture if the product are storde at 2-8°C.

Ref. No.	description
ZR106020	20 x 900 μl solution