

THE NEW CRYOTOP®

NOT ONLY THE BEST PROVEN RESULTS BUT ALSO THE BEST USER EXPERIENCE

KITAZATO®

NEW VERSION



OLD VERSION




WIDE STRIP TO LOAD THE SPECIMEN

The strip's been widened. Now, with a bigger surface, it is even easier to load the specimens.



EASIER INSERTION

The diamond tip and the reshaped conjunction between strip and handle grant a smoother insertion into the cover straw.



IDENTIFICATION MARKS

The Cryotop® has now two identification marks on its upper surface, defining the side where the specimens should be loaded, if you are vitrifying, or are loaded, if you are thawing.



DIBIMED

C/Luis Buñuel, 1 Oficinas Pta. 4
46015 Valencia Spain



www.kitazato-dibimed.com