

## Certificate of Analysis GM508 CultActive

Batchnumber	CACT-1803
Production date	19 April 2018
Expiry date	31 October 2018
Total volume (mL)	990

	Results	Specifications
pH	7,29	7,00-7,60
Osmolality	420	400-450 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,025	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	95,24	≥ 80 %

Release date 04 May 2018

Released by Dr. Julia Heinzmann

Signature

*p.p. Julia Heinzmann*