

**Certificate of Analysis
GM501 Cult with Gentamicin**

Batchnumber	CULT-1806G
Production date	16 August 2018
Expiry date	28 February 2019
Total volume (mL)	33100
Batchnumber HSA	K809B6671

	Results	Specifications
pH	7,26	7,20-7,45
Osmolality	278	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,05	< 0,25 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

Release date 03 September 2018

Released by Dr. Julia Heinzmann

Signature

p.p. Julia Heinzmann