

Certificate of Analysis

Product name: OOSAFE® HAND DISINFECTANT

Product code: OOHD

Batch/LOT No.: 81010

Production date: 2021-11

Expiry date: 2024-11

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid	Clear liquid	NA
Odor / Olfactive	Mild	Mild	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	8,0 - 9,5	9.2	NA
Density, 20 °C, gr/cm ³ / Density measuring device	1,01 – 1,10	1,08	NA
Quaternary Ammonium Compounds / QUAT test stripe	2,4 – 2,8	2,6 mg / gr	NA
Additives		8,2 mg / gr	NA
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-7514-0122	Passed	31/01/2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5611-0122SPAR	Passed	28/01/2022

Quality Control according to ISO 13458:2016


GOosafe with SparMED!


Andrius Daugintis

Quality Control Department

SparMED ApS

Date: 01/02/2022



 <p>GOosafe in IVF Laboratories</p>	PREPARED BY:	APPROVED BY:
	QUALITY CONTROL MANAGER	MANAGEMENT REPRESENTATIVE
	Signature <i>lvei</i>	Signature <i>Bmfe</i>