

## Certificate of Analysis

**Product name: Oosafe® Hand Disinfectant****Product code: OOHD****Batch/LOT No.: 81030****Production date: 2022-11****Expiry date: 2025-11**

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid	Clear liquid	N/A
Odor / Olfactive	mild	mild	N/A
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	8,0 – 9,5	9,4	N/A
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1,01 – 1,10	1,08	N/A
Quaternary Ammonium Compounds / QUAT test stripe	2,4 – 2,8	2,6 mg/gr	N/A
Additives		8,2 mg/gr	N/A
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR – 1937 – 0523	Passed	23/05/2023
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S6527 – 0523SPAR	Passed	26/05/2023

Quality Control according to ISO 13485:2016

**GOosafe with SparMED!**

Lana Shenkevych

Quality Control Department

SparMED ApS

Date: 27/05/2023

