

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber OIL-1834
 Production date 08 November 2018
 Expiry date 31 May 2020
 Total volume (ml) 183000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	21,2	< 30 cP at 30°C
Specific gravity	0,8431	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	90,48	≥ 80 %

Release date 26 November 2018

Released by Vanessa Bierlich

Signature