

**Certificate of Analysis  
GM501 Mineral Oil**

Batchnumber	OIL-1838
Production date	11 December 2018
Expiry date	30 June 2020
Total volume (ml)	199000

	<b>Results</b>	<b>Specifications</b>
Appearance	PASSED	PASSED
Sterility test <small>Sterilised by sterile filtration over 0,2µm filter</small>	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	22,3	< 30 cP at 30°C
Specific gravity	0,8516	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay <small>(1-cell to expanded blastocyst after 120 hours)</small>	95,24	≥ 80 %

Release date 02 January 2019

Released by Vanessa Bierlich

Signature