

**Certificate of Analysis
GM501 Bromelain in Dulbecco's PBS**

Batchnumber	BROM-2003
Production date	01 December 2020
Expiry date	30 June 2021
Total volume (mL)	2200

	Results	Specifications
pH	7,32	7,20-7,50
Osmolality	282	270-290 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Endotoxin	< 0,097	< 0,50 EU/mL
Mouse Embryo Assay (1-cell)	100	≥ 80 %
Sperm Survival Test	96	≥ 80 %

Release date	16 December 2020
Released by	Vanessa Bierlich

Signature