

RECORDING MODUL

CAMERAS: 1MOS HD / 3MOS HD / 4K UHD



Technical Data

- Video format:
- Full HD 1920 x 1080 p
 - Mpeg4 (H.264 Codec)
 - Frame rate: 60fps
- Image format:
- JPG (1920 x 1080)

Recording time (approximately)

H.264	1920x1080p / 60fps	1920x1080p / 60fps	1920x1080p / 60fps
BIT Rate	24 Mbps (Best quality)	12 Mbps	4Mbps
USB Stick - 4GB	22 Minutes	45 Minutes	2 Hours 15 Minutes
USB Stick - 16GB	1 Hour 30 Minutes	3 Hours	9 Hours
USB Stick - 32GB	3 Hours	6 Hours	18 Hours
USB Stick - 64GB	6 Hours	12 Hours	36 Hours
USB Stick - 128GB	12 Hours	24 Hours	72 Hours
USB Stick - 256GB	24 Hours	48 Hours	144 Hours

