


Certificate of Analysis GM501 GentleVit Thaw

Batchnumber	GT-2506
Production date	27 February 2025
Expiry date	28 February 2026
Total batch size	201
Batchnumber HSA	B03H108433

Parameter	Results	Specifications
GM501 GentleVit Thaw 1		
pH	7,26	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Concentration HSA	16,3	13,9-20,9 g/l
GM501 GentleVit Thaw 2		
pH	7,29	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Concentration HSA	16,9	14,4-21,6 g/l
GM501 GentleVit Thaw 3		
pH	7,26	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	824	805-850 mOsm/kg
Concentration HSA	17,1	14,9-22,3 g/l
GM501 GentleVit Thaw 4		
pH	7,29	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	545	535-565 mOsm/kg
Concentration HSA	18,4	15,4-23,1 g/l
GM501 GentleVit Thaw 5		
pH	7,29	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	414	405-435 mOsm/kg
Concentration HSA	19,3	15,7-23,6 g/l
GM501 GentleVit Thaw 6		
pH	7,29	7,20-7,40
Sterility test ¹	No growth	No growth
Endotoxin	< 0,1	< 0,25 EU/mL
Osmolality	279	270-290 mOsm/kg
Concentration HSA	19,9	16,0-24,0 g/l
Mouse Embryo Assay (1-cell)	100	≥ 80 %

Source materials from which this product was derived was found negative when tested for antibodies against HIV-1/2 and HCV and non-reactive for HbsAg, HBV RNA, HIV-1 RNA and syphilis. No known test method can offer assurance that products derived from human blood will not transmit infectious agents. Therefore, handle all specimens as if capable of transmitting HIV or hepatitis.

¹Sterility test (according to the current PH. Eur. 2.6.1 / USP <71>)

Release date	1. April 2025
Released by	Nicole Geiger
Signature	

REV 02_00 / 2023-09-15