

Certificate of Analysis
Denu-Tips

Batchnumber 44190103
Expiry date 31 January 2022

Article No.

GY-135/xx
GY-140/xx
GY-145/xx
GY-150/xx
GY-165/xx
GY-170/xx
GY-175/xx
GY-200/xx
GY-275/xx
GY-300/xx

	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	95,24	≥ 80 %

Release date 11 March 2019
Released by Dr. Julia Heinzmann
Signature

p.p. Julia Heinzmann