

**Certificate of Analysis**  
**Denu-Tips**

Batchnumber 442303/34/01  
 Production date 21 August 2023  
 Expiry date 31 July 2026

**Article No.**

GY-170/20

<b>Parameter</b>	<b>Result</b>	<b>Specification</b>
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	95	≥ 80 %

Release date 02 October 2023  
 Released by Nicole Geiger  
 Signature

pp. N. Geiger

REV 08\_00 / 2023-09-15