

**Certificate of Analysis  
GM501 Collagenase**

Batchnumber	COLL-1904
Production date	10 July 2019
Expiry date	30 April 2020
Total volume (mL)	1130

	<b>Results</b>	<b>Specifications</b>
pH	7,43	7,10-7,50
Osmolality	280	270-310 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Sperm survival Test 4 hours exposure	88	≥ 80 %

Release date	25.07.219
Released by	Dr. Julia Heinzmann
Signature	

*p.p. Julia Heinzmann*