

**Certificate of Analysis
GM501 Collagenase**

Batchnumber	COLL-1906
Production date	22 November 2019
Expiry date	31 August 2020
Total volume (mL)	910

	Results	Specifications
pH	7,34	7,10-7,50
Osmolality	285	270-310 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Sperm survival Test 4 hours exposure	95	≥ 80 %

Release date	10 December 2019
Released by	Vanessa Bierlich
Signature	

P.P. 

