

**Certificate of Analysis  
GM501 Collagenase**

Batchnumber	COLL-2005
Production date	03 September 2020
Expiry date	30 September 2021
Total volume (mL)	1000

	<b>Results</b>	<b>Specifications</b>
pH	7,35	7,20-7,50
Osmolality	287	270-310 mOsm/kg
Bacterial screen	PASSED	PASSED
Fungal screen	PASSED	PASSED
Sperm survival Test 4 hours exposure	93	≥ 80 %

Release date	21 September 2020
Released by	Vanessa Bierlich
Signature	

P.P. 