

**Certificate of Analysis
GM501 Collagenase**

Batchnumber COLL-2402
 Production date 02 May 2024
 Expiry date 28 February 2025
 Total volume (ml) 1330

Parameter	Results	Specifications
pH	7,43	7,20-7,70
Osmolality	281	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	NO GROWTH	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	97	≥ 80 %

Release date 17 May 2024
 Released by Nicole Geiger
 Signature

