

**Certificate of Analysis
GM501 Collagenase**

Batchnumber COLL-2403
 Production date 02 July 2024
 Expiry date 30 April 2025
 Total volume (ml) 1670

Parameter	Results	Specifications
pH	7,41	7,20-7,70
Osmolality	283	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	91	≥ 80 %

Release date 17 July 2024
 Released by Nicole Geiger
 Signature
