

## Certificate of Analysis GM501 Collagenase

Batchnumber	COLL-2404
Production date	13 September 2024
Expiry date	30 June 2025
Total volume (ml)	2110

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,31	7,20-7,70
Osmolality	285	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	NO GROWTH	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	92	≥ 80 %

Release date	30 September 2024
Released by	Nicole Geiger
Signature	