

**Certificate of Analysis
GM501 Collagenase**

Batchnumber COLL-2405
 Production date 03 December 2024
 Expiry date 30 September 2025
 Total volume (ml) 2790

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,4	7,20-7,70
Osmolality	279	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	91	≥ 80 %

Release date 18 December 2024
 Released by Nicole Geiger
 Signature

