

**Certificate of Analysis
GM501 Collagenase**

Batchnumber COLL-2601
 Production date 02 March 2026
 Expiry date 31 December 2026
 Total volume (ml) 1560

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,27	7,20-7,70
Osmolality	279	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	95	≥ 80 %

Release date 17 March 2026
 Released by Dr. René Richarz
 Signature

