

Certificate of Analysis GM501 Collagenase

Batchnumber	COLL-2602
Production date	04 May 2026
Expiry date	28 February 2027
Total volume (ml)	1560

Product was sterilized by aseptic processing techniques

Parameter	Results	Specifications
pH	7,34	7,20-7,70
Osmolality	279	270-310 mOsm/kg
Sterility test (according to the current Ph.Eur.2.6.1 USP <71>)	No growth	No growth
HSSA (% motility compared with control after 4 h exposure to test medium)	100	≥ 80 %

Release date	20 May 2026
Released by	Nicole Geiger
Signature	

REV 11_00 / 2023-09-15