

KITAZATO®

Now you have more
time to *freeze time*



The new VT 600

*Vitrification and Thawing media for the Cryotop® Method
Now with 12 months shelf life.*

In Kitzato we operate under the **Constant Improvement Philosophy**. The constant feedback from IVF professionals is essential to us.

A brand that listens to its clients is a brand that improves its methods.

We are pleased to launch the **New VT600 series**, the media that answers your storage needs with a **shelf life of 12 months**.

So, no need to rush, you can keep those vitrification and thawing procedures coming!

Evidence of success

The Cryotop® Method media has evolved to continue offering outstanding results.

Our Method is backed up by numerous independent clinical studies with the biggest sample

Oocytes:



n=3286
Cobo
2010



n=124
Rienzi
2010

Blastocysts:



n=2543
Shi
2012



n=6861
Liu
2013

Embryos:



n=2031
Cobo
2012



n=553
Rodríguez-Purata
2010



Choose Science

Kitazato always lead vitrification with innovation. This is what makes our media special.

- **We were the first ones to replace sucrose with trehalose**, demonstrating its superiority as an osmotic agent.
- **Our media has evolved replacing the HSA (Human Serum Albumin) for HPC (Hydroxypropyl Cellulose)**, reducing the mechanical stress during thawing.
- **The incorporation of Gentamicin** prolongs the sterility of the media, guaranteeing greater safety in handling the solutions.
- **The HEPES buffer and the sodium bicarbonate act as stabilizers for the pH**, ensuring normal cellular physiological conditions during vitrification and thawing procedures.