

**Certificate of Analysis
GM501 Mineral Oil**

Batchnumber	OIL-2018
Production date	13 July 2020
Expiry date	31 January 2022
Total volume (mL)	180000

	Results	Specifications
Appearance	PASSED	PASSED
Sterility test Sterilised by sterile filtration over 0,2µm filter	PASSED	PASSED
Acidity	Conform Ph Eur	Conform Ph Eur
Solid paraffin	Conform Ph Eur	Conform Ph Eur
Viscosity	24,6	< 30 cP at 30°C
Specific gravity	0,8429	< 0,86 g/ml
Endotoxin	< 0,1	< 0,1 EU/ml
Peroxide value (POV)	< 0,1	< 0,1 mEq/kg
Mouse Embryo Assay (1-cell to expanded blastocyst after 120 hours)	100	≥ 80 %

Release date 03 August 2020

Released by Ileana V. Hinz

Signature

