

**Certificate of Analysis  
Micropipettes - Biopsy**

Batchnumber 20220926  
 Production date 30 September 2022  
 Expiry date 30 September 2027

**Article No.**  
 BP-20-35°-bl

	<b>Results</b>	<b>Specifications</b>
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	90	≥ 80 %

Release date 05 January 2023  
 Released by Nicole Geiger

Signature

