

**Certificate of Analysis
Micropipettes - Holding**

Batchnumber 20221215
 Production date 15 December 2022
 Expiry date 31 December 2027

Article No.

HP-100-30°
 HP-100-35°

	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	95	≥ 80 %

Release date 31 January 2023

Released by Nicole Geiger

Signature



REV 03_00 / 2020-01-27

