

**Certificate of Analysis
Micropipettes - ICSI**


Batchnumber 20240220
 Production date 28 February 2024
 Expiry date 28 February 2029

Article No.

IC-5-35°-s
 IC-SL-4.6-30°-s

Parameter	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	100	≥ 80 %

Release date 17 April 2024
 Released by Nicole Geiger
 Signature



REV 05_00 / 2023-09-15