

**Certificate of Analysis
Micropipettes - Holding**

Batchnumber 20250131
 Production date 31 January 2025
 Expiry date 31 January 2023

Article No.


HP-120-30°

Parameter	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	100	≥ 80 %

Release date 04 March 2025

Released by Nicole Geiger

Signature



REV 04_00 / 2023-09-15