

Certificate of Analysis Micropipettes - Biopsy

Batchnumber 20250220

Expiry date 28 February 2030

Article No.

BP-30-30°-bl

BP-35-30°-bl

Parameter	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	100	≥ 80 %

Release date 31 March 2025

Released by Nicole Geiger

Signature

pp. N. Geiger

REV 04_00 / 2023-09-15