


Certificate of Analysis Micropipettes - Holding

Batchnumber	20251007
Production date	07 October 2025
Expiry date	30 October 2030

Article No.

HP-80-30°
HP-100-30°

Parameter	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	100	≥ 80 %

Release date	13 November 2025
Released by	Nicole Geiger
Signature	

REV 04_00 / 2023-09-15