

Certificate of Analysis Micropipettes - Holding

Batchnumber		20260409
Production date		09 April 2026
Expiry date		25 March 2031

Article No.

HP-100-30°
HP-100-35°
HP-120-30°

Parameter	Results	Specifications
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	100	≥ 80 %

Release date		01 June 2026
Released by		Nicole Geiger
Signature		

pp. N. Geiger

REV 04_00 / 2023-09-15