

Certificate of Analysis

Product name: Oosafe® Disinfection Wipes (Medical Devise)

Product code: OOMDW-120

Batch/LOT No.: 21002-MD

Production date:2022-03

Expiry date:2025-03

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid inside wipes	Clear liquid inside wipes	NA
Odor / Olfactive	Odorless	Odorless	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	7.0 – 9.0	8.3	NA
Density, 20 °C, gr/cm ³ / Density measuring device	1.001 – 1.010	1.008	NA
Quaternary Ammonium Compounds / QUAT test stripe	1,9 – 2,2 %	2.0 %	NA
Water		the rest	NA
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-8096-0422	Passed	23-04-2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5779-0422SPAR	Passed	20-04-2022

Quality Control according to ISO 13458:2016

GOosafe with SparMED!

Zahraa Baker
Quality Control Department
SparMED ApS
Date: 09-05-2022

