

# Certificate of Analysis

**Product name: OOSAFE® HAND DISINFECTANT**

**Product code: OODH**

**Batch/LOT No.: 81010**

**Production date: 2021-11**

**Expiry date: 2024-11**

Analysis / Method	Specification / Reference No.	Result	Completion date
Appearance / Visual	Clear liquid	Clear liquid	NA
Odor / Olfactive	Mild	Mild	NA
pH-Rate (2% active substance in water) / pH meter (DIN EN 1262)	8,0 - 9,5	9.2	NA
Density, 20 °C, gr/cm <sup>3</sup> / Density measuring device	1,01 – 1,10	1,08	NA
Quaternary Ammonium Compounds / QUAT test stripe	2,4 – 2,8	2,6 mg / gr	NA
Additives		8,2 mg / gr	NA
Water		The rest	
Proven non-embryotoxic by Mouse Embryo Assay Test. Over 80% embryo development to the expanded blastocyst stage within 96 hours.	SPAR-7753-0222	Passed	07/03/2022
Proved stable human sperm motility: ≥75% sperm motility after 24 hours proven.	S5690-0322SPAR	Passed	03/03/2022

Quality Control according to ISO 13458:2016


**GOosafe with SparMED!**


Andrius Daugintis

Quality Control Department

SparMED ApS

**Date: 08/03/2022**



 <p>GOosafe in IVF Laboratories</p>	PREPARED BY:	APPROVED BY:
	QUALITY CONTROL MANAGER	MANAGEMENT REPRESENTATIVE
	Signature <i>lvei</i>	Signature <i>Bmbe</i>