



Certificate of Analysis

Product: SpermFilter® Gradient 45%

Batch number: FLT45-209
Expiration Date: 30-09-2023
(DD/MM/YYYY)

Analysis	Limits	Results
pH:	7,20 – 7,60	7,35
Osmolality:	310 – 340 mOsm/kg	323
Sterility Test Bacteria:	NEGATIVE	PASSED
Sterility Test Fungi:	NEGATIVE	PASSED
(Sterile by filtration over 0,2µm filter)		
Endotoxin:	< 0,5 EU/ml	< 0,100 EU/ml
Sperm Survival:	≥ 80%	99%
(4 hours exposure)		
Sperm Survival:	≥ 75%	94%
(24 hours exposure)		

Approved on: 07-04-2022
(DD/MM/YYYY)
By: M. Berbers

Signature:

Quality Control Department, Gynotec B.V.