



Kitazato Corporation


<http://www.kitazato.co.jp/>

PRODUCT: Hypure® Oil Heavy

LOT#: D0201

REF: 96011  
Code: MOHVD-100P

EXPIRY DATE : 2024/08/31  
Year/Month/Day

  
APPROVED BY : Kazumi Shibata  
Quality Control  
2023/02/16

### **CERTIFICATE OF ANALYSIS**

The product has been manufactured by Kitazato Corporation in conformity with requirements of ISO13485 Quality management system. The quality of products was tested and evaluated by quality control procedures

ASSAY	SPECIFICATION	RESULT
Color:	Clear, Particle Free	PASSED
Endotoxin:	< 0,1 EU/ml	<0.1 EU/mL
Sterility:	No microbial Growth	PASSED
MEA <sup>1</sup> :		
- Control:	≥ 80%	100%
- This Lot:	≥ 80%	90.48%
POV:	<0,1 mEq/Kg	<0,1 mEq/kg

<sup>1</sup> One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing to fully expanded blastocysts after 120 hours in culture medium.