



Kitazato Corporation

<http://www.kitazato.co.jp/>

PRODUCT: Hypure™ Oil Heavy

LOT#: B0921

REF: 96011
Code: MOHVD-100P

EXPIRY DATE : 2023/03/31

Year/Month/Day

APPROVED BY : Kazumi Shibata
Quality Assurance

2021/10/12

CERTIFICATE OF ANALYSIS

The product has been manufactured by Kitazato Corporation in conformity with requirements of ISO13485 Quality management system. The quality of products was tested and evaluated by quality control procedures

ASSAY	SPECIFICATION	RESULT
Color:	Clear, Particle Free	PASSED
Endotoxin:	< 0,1 EU/ml	<0.100
Sterility:	No microbial Growth	PASSED
MEA ¹ :		
- Control:	≥ 80%	93,33%
- This Lot:	≥ 80%	95,24%
POV:	<0,1 mEq/Kg	<0.1

¹ One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing to fully expanded blastocysts after 120 hours in culture medium.