

PRODUCT: Hypure™ Oil Heavy

LOT#: C0511

REF: 96015
Code: MOHVD-500P

EXPIRY DATE : 2023/11/30

Year/Month/Day



APPROVED BY : Kazumi Shibata
Quality Assurance

2022/05/27

CERTIFICATE OF ANALYSIS

The product has been manufactured by Kitazato Corporation in conformity with requirements of ISO13485 Quality management system. The quality of products was tested and evaluated by quality control procedures

| ASSAY | SPECIFICATION | RESULT |
|--------------------|----------------------|--------|
| Color: | Clear, Particle Free | PASSED |
| Endotoxin: | < 0,1 EU/ml | <0.100 |
| Sterility: | No microbial Growth | PASSED |
| MEA ¹ : | | |
| - Control: | ≥ 80% | 86,67% |
| - This Lot: | ≥ 80% | 95,24% |
| POV: | <0,1 mEq/Kg | <0.1 |

¹ One cell mouse embryo testing is performed using fresh mouse embryos. Test results indicate the percentage of mouse embryos developing to fully expanded blastocysts after 120 hours in culture medium.