


PREMIUM FAS SINGLE LUMEN

FOLLICLE ASPIRATION SET

The Premium FAS, the follicle aspiration set is implemented for ultrasound-guided, transvaginal harvesting of oocytes from follicles for IVF and ICSI.

FEATURES


- standard follicle aspiration needle with regular tubing
- sharp lancet and back cut bevel tip available that ensure a clean facet cutting
- diamond laser ultrasound markings on tip which offer stable and excellent vision during ultrasound
- translucent tubing in the same size as needle for safe oocyte retrieval
- arrow marking on handle for bevel point indication
- handle shape facilitates rotation



ORIGINAL PICTURE

GYNO INFO PREMIUM FAS SINGLE LUMEN

Sterile	sterilized by gamma irradiation	
Packaging	50/shipping box (standard configuration)	
Quality control	MEA and LAL	
Order number	#4551 NS-AS1	Ⓜ
	#4551 NS-BS1	Ⓜ
	#4551 NS-CS1	Ⓜ
	#4551 NS-DS1	Ⓜ
	#4551 NS-AS1BCB	Ⓜ
	#4551 NS-CS1BCB	Ⓜ
	#4551 NS-DS1BCB	Ⓜ



Gynetics Medical Products n.v.
 Rembert Dodoensstraat 51 / 3920 Lommel, Belgium
 tel. +32 11 64 58 72 / fax +32 11 64 06 77
 info@gynetics.be / www.gynetics.be

GYNO INFO

PREMIUM FAS SINGLE LUMEN

① Order number	Gauge	OD/ID (mm)	Length (mm)	Working length (mm)	Bevel	Control
#4551 NS-AS1	17G	1.4/1.2	400	350	Lancet	Flexible
#4551 NS-BS1	17G	1.4/1.2	340	290	Lancet	Flexible
#4551 NS-CS1	17G	1.4/1.0	400	350	Lancet	Rigid
#4551 NS-DS1	17G	1.4/1.0	340	290	Lancet	Rigid
#4551 NS-AS1BCB	17G	1.4/1.2	400	350	Back cut bevel	Flexible
#4551 NS-CS1BCB	17G	1.4/1.0	400	350	Back cut bevel	Rigid
#4551 NS-DS1BCB	17G	1.4/1.0	340	290	Back cut bevel	Rigid

Order number	Gauge	OD/ID (inch)	Length (inch)	Working length (inch)	Bevel	Control
#4551 NS-AS1	17G	0.055/0.047	15.75	13.779	Lancet	Flexible
#4551 NS-BS1	17G	0.055/0.047	13.39	11.417	Lancet	Flexible
#4551 NS-CS1	17G	0.055/0.039	15.75	13.779	Lancet	Rigid
#4551 NS-DS1	17G	0.055/0.039	13.39	11.417	Lancet	Rigid
#4551 NS-AS1BCB	17G	0.055/0.047	15.75	13.779	Back cut bevel	Flexible
#4551 NS-CS1BCB	17G	0.055/0.039	15.75	13.779	Back cut bevel	Rigid
#4551 NS-DS1BCB	17G	0.055/0.039	13.39	11.417	Back cut bevel	Rigid

- ② Tubing from needle to cork : Regular tubing 100 cm / 39.369 inch
- ③ Size silicone cork : Small 14 x 18 x 20 mm / 0.551 x 0.709 x 0.787 inch
- ④ Tubing from cork to pump : Regular tubing 100 cm / 39.369 inch

