

Certificate of Analysis Oocyte Aspiration Needle

Batchnumber	712603/33/01
Production date	10 August 2026
Expiry date	31 July 2029

Article No.

GY-711724LF

Parameter	Result	Specification
Packaging	PASSED	PASSED
Sterilization	PASSED	PASSED
Mouse Embryo Assay (1-cell, blastocysts after 96 hours)	86	≥ 80 %

Release date	29 October 2026
Released by	Dr. René Richarz
Signature	



REV 04_00 / 2025-05-21