

GYNEMED

Helps when others fail



www.gynemed.de



Reproductive medical solutions by Gynemed

- Oocyte Handling
- Culture
- Sperm Processing
- Cryopreservation
- In-Vitro-Diagnostics
- Laboratory Equipment
- Diagnostics
- Consumables



CultActive

GM508 CultActive is a bicarbonate buffered HSA-free reagent designed to investigate oocytes of patients with failed fertilization after previous intracytoplasmic sperm injection cycles. GM508 CultActive is designed to investigate if fertilization failure after previous ICSI cycles is due to a deficient oocyte activation.



Collagenase

GM501 Collagenase is a reagent for the digestion of human testicular tissue obtained by testicular biopsy (TESE) for in-vitro examination procedures.

Using Collagenase, it is possible to degrade testicular tissue in single cells to facilitate the isolation of sperm cells.



SpermMobil

GM501 SpermMobil is a HEPES buffered HSA-free reagent containing low bicarbonate.

It is used for in-vitro examination of sperm cells of necrozoospermic ejaculates as well as of immotile sperms isolated from testicular tissue (TESE).



Bromelain in Dulbecco's PBS

Bromelain in Dulbecco's PBS is designed for liquefaction of viscous semen samples prior to semen analysis and preparation for further IVF treatment.



GM501 MediaLine

From routine to special application

Oocyte Handling

Culture

Sperm Processing

Cryopreservation

In-Vitro-Diagnostics



DENU-Tips

Made out of Polyamid (BPA-free)

Vials of 20

Sterilized for 3 years

Different sizes available (80 - 550 μm)


Mouse-Embryo-tested

CE marked



Pipettes for Micromanipulation

Made out of borosilicate glass
Customized production available
Mouse-Embryo-tested
CE marked



Oocyte Aspiration Needles


Single or double lumen
Echogenic markings
Handles with position marking
Variable needle and tubing configurations
Mouse-Embryo-tested
CE marked

IVFtech



Laboratory Equipment

We additionally offer laboratory devices and consumables by renowned brands such as IVFtech, Hamilton Thorne, Planer and others as exclusive distributor including service and maintenance.



TrakJector™

All-in-one Trackball or Joystick offers true one-handed operation.

Available in two injector visions:

MIU^{mot} Oil or MIU^{mot} AIR

Fits on all common inverted microscopes

Gynemed has been dedicated to the development, manufacture and distribution of high-quality medical products for reproductive medicine since 1998.

Daily contact to our customers enables us to support them in their goal to offer the best medical care possible in the field of assisted reproduction.

GYNEMED

a Hamilton Thorne Ltd. company

Gynemed GmbH & Co. KG

Lübecker Straße 9
23738 Lensahn

Germany

Tel.: +49 (0) 4363 90 32 90

Fax: +49 (0) 4363 90 32 9-19

info@gynemed.de

www.gynemed.de

