


# GYNEMED

## Assisted Reproduction



Oocyte Aspiration

# Oocyte Aspiration Needles

## KEY features

- **Sterile** - Sterilized by gamma irradiation and individually packed.
- **Echogenic markings** - Enable excellent visualization under ultrasound view.
- **Customized production** - Available upon request.
- **CE marked** - To meet international standards (Class IIa).
- **Shelf life** - Sterility and functionality guaranteed for 3 years from time of manufacture.
- **MEA-tested** - A certificate of analysis (COA) is available for each lot number upon request from our website.

## Oocyte Aspiration Needle

The oocyte aspiration needles are available as single or double lumen sets.





The single lumen needles are available with or without flushing (3-way-valve).

The double lumen needles can be used for concurrent aspiration and flushing (stopcock directly connected to handle).

The silicone stopper is suitable for 14 ml-round bottom collection tubes.

The oocyte aspiration needles have comfortable handles with position marking.

Variable needle and tubing configurations are available (length, thickness etc.).



# Oocyte Aspiration Needle Specifications

REF	Product description	Unit
<b>GY-711724LF</b>	single lumen, without flushing, 17 gauge needle, length of the needle 32 cm, back cut bevel, echomark microblasting, length of aspiration tubing 100 cm, length of vacuum tubing 100 cm, luer female	10
<b>GY-711724LF-A</b>	single lumen, without flushing, 17 gauge needle, length of the needle 35 cm, back cut bevel, echomark microblasting, length of aspiration tubing 100 cm, length of vacuum tubing 100 cm, luer female	10
<b>GY-711824LF</b>	single lumen, without flushing, 18 gauge needle, length of the needle 32 cm, back cut bevel, echomark microblasting, length of aspiration tubing 100 cm, length of vacuum tubing 100 cm, luer female	10
<b>GY-721724LF</b>	single lumen, with flushing, 17 gauge needle, length of the needle 32 cm, back cut bevel, echomark microblasting, length of aspiration tubing 100 cm, length of vacuum tubing 100 cm, luer female	10
<b>GY-931627LF-A</b>	double lumen, with flushing, outer diameter 16 gauge needle, length of the needle 35 cm, back cut bevel, echomark dots up to the end of the tips, length of aspiration tubing 100 cm, length of vacuum tubing 100 cm, luer female	10

*Other specifications available upon request.*

Designed and distributed by:



Gynemed  
GmbH & Co. KG  
Lübecker Straße 9  
23738 Lensahn  
Germany

Tel.: +49 (0) 4363 90 32 9-0  
Fax: +49 (0) 4363 90 32 9-19

info@gynemed.de  
www.gynemed.de

Hamilton Thorne, Inc.  
100 Cummings Center  
Suite 465E  
Beverly, MA 01915  
USA

Tel.: +1 978-921-2050  
Tel.: +1 800-323-0503  
Fax: +1 978-921-0250

sales@hamiltonthorne.com  
www.hamiltonthorne.com

- Made in Germany -