

GYNEMED



Years



WITH **GYNÉTICS**[®] A NEW GYNEMED BRAND!

www.gynemed.de

GM501 MediaLine

From routine to special application

Oocyte Handling

Culture

Sperm Processing

Cryopreservation

In-Vitro-Diagnostics



Unique solutions



CultActive

GM508 CultActive is a bicarbonate buffered HSA-free reagent designed to investigate oocytes of patients with failed fertilization after previous intracytoplasmic sperm injection cycles. GM508 CultActive is designed to investigate if fertilization failure after previous ICSI cycles is due to deficient oocyte activation.



Collagenase

GM501 Collagenase is a reagent for the digestion of human testicular tissue obtained by testicular biopsy (TESE).

Using Collagenase, it is possible to degrade testicular tissue in to single cells to facilitate the isolation of sperm cells for use in ART procednes.



SpermMobil

GM501 SpermMobil is a HEPES buffered HSA-free reagent containing low bicarbonate.

It is used for in-vitro examination of sperm cells of necrozoospermic ejaculates as well as of immotile sperms isolated from testicular tissue (TESE).



Bromelain in Dulbecco's PBS

Bromelain in Dulbecco's PBS is designed for liquefaction of viscous semen samples prior to semen analysis and preparation for further IVF treatment.



Oocyte Aspiration Needles

Single or double lumen

Echogenic markings

Handle shape facilitates rotation

Variable needle and tubing configurations

Mouse-Embryo-Assay

CE marked



Catheters

Embryo Transfer Catheters
for routine, trial and difficult transfer

Insemination Catheters
with or without memo function


Mouse-Embryo-Assay

CE marked




Pipettes for Micromanipulation

Made out of borosilicate glass
Customized production available
Mouse-Embryo-Assay
CE marked



Strip-ID

Flexible and unbreakable
material that will not scratch dishes
Vials of 10
Different sizes available (125 - 300 µm)
Mouse-Embryo-Assay
CE marked



GyPetto & GyPettino

Gynemed Pipettor

Easy disinfection

No lacquers or colours applied

Two sizes available

IVFtech



Laboratory Equipment

We additionally offer laboratory devices and consumables by renowned brands such as IVFTech, Hamilton Thorne, Planer and others as exclusive distributor including service and maintenance.

Gynemed has been dedicated to the development, manufacture and distribution of high-quality medical products for reproductive medicine since 1998.

Daily contact to our customers enables us to support them in their goal to offer the best medical care possible in the field of assisted reproduction.

GYNEMED

a Hamilton Thorne Ltd. company

info@gynemed.de
www.gynemed.de

