

# CERTIFICATE OF ANALYSIS

<b>Product</b>	100 µm chamber depth - 2 analysis chambers	
<b>Slide reference</b>	SC 100-01-02-A-CE	
<b>Batch   Machine</b>	042553	L
<b>Production date range (DD-MM-YYYY)</b>	20-01-2025 - 20-01-2025	
<b>Best before (Month-YY)</b>	January-26	

<b>Analysis</b>	<b>Criteria</b>	<b>Result</b>
Dimensions and microscope inspection	20 inspection points	PASSED
Semen Toxicity Test (internal protocol): <i>Hamilton Thorne IVOS II Digital</i> Type of semen: <i>Porcine</i>	Time degradation 4 minutes <8%	PASSED
Chamber height: <i>Low Coherence Interferometry</i>	Median <10% from indicated mean. Variation <2%	PASSED
Average chamber height (µm)		99,40

**WE CERTIFY THAT THE ABOVE MENTIONED PRODUCTS  
ARE IN COMPLIANCE WITH OUR SPECIFICATIONS**

<b>Date</b>	January 20, 2025
<b>Name</b>	Anne Linda van Kappel
<b>Qualification</b>	General Manager
<b>Signature</b>	