

CERTIFICATE OF ANALYSIS

Product	accu-beads®+	Part Number	710257
----------------	--------------	--------------------	--------

Concentration Type	Acceptable Range	Lot #		Date of Manufacture	Expiration Date
18 M/ml	± 2.5 M/ml	241211181		2024.11.12	2026.11.13
Test Type and Quantity	20 µm Fixed Coverslip Chamber Calibration	S.D.	CV%	Calculated Hemacytometer Calibration Range*	
1 L - Bio	18.67	0.23	1.24	23.97-24.57	
1 L - QC	18.79	0.39	2.06	23.29-24.93	
5 ml - QC	18.76	0.19	1.02	24.14-24.64	
QC Analysis by: <i>Carlos A. Ulloa</i>		Reviewed by: Nicole Thibodeau			
Date: <i>2024-11-22</i>		Date: 2024-11-22			
Signature: <i>Carlos A. Ulloa</i>		Signature: Nicole Thibodeau <small>Digitally signed by Nicole Thibodeau Date: 2024.11.22 16:30:51 -05'00'</small>			

* Calculation based on hemacytometer flow dynamics ratio as described in accu-beads Research Report 45. Please contact Hamilton Thorne for more information.